ADNucleis

ARCHITECT IN DIAGNOSIS

CATALOGUE

MATERIALS ET CONSOMABLES

February 2016
V7

ADNUCLEIS SAS

1 246 000 € Capital

Head Offices : 30 chemin des Mouilles
69290 GREZIEU LA VARENNE - FRANCE

SIRET : 497 802 165 000 17

TVA intraCEE : FR 3649780216500017

Email : contact@adnucleis.com

Tel : + 33 (0) 478 571 640

Mob : + 33 (0) 660 960 141

TABLE of CONTENTS

I.	Conditions of a PCR lab implementation	3
1	Necessary material.....	3
2	Training.....	5
3	Prix des kits ADNucleis	6
4	Avantages des kits ADNucleis	6
II.	Methodologie.....	6
1	Enrichissement et milieux de culture.....	6
2	Extraction et purification d'ADN	7
3	Les kits d'amplification.....	7
4	External Point Control (EPC) & Internal Point Control (IPC)	8
III.	Products.....	9

I. Conditions of a PCR lab implementation

1 Necessary material

Real time RT PCR or QPCR

APPLIED BIOSYSTEMS
7000-7300-7500-7900

**STRATAGENE-
AGILENT**
MX 3005P

QIAGEN
Rotorgene

EPPENDORF
Realplex

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THERMOCYCLER

**MULTICHANNEL
MICROPIPETTE**
1-10 μ L / 30-50 μ L

**PLATE OR STRIPS
CENTRIFUGE**

STOMACHER

**MICROBIOLOGICAL SAFETY CABINET
(MSC-P2)**

REFRIGERATOR / FREEZER

Material installation + trials : 48h.

Level P1 laboratory : no specific sanitary requirements.

Ventilation : air extraction from the amplification room /air entry inside the DNA extraction room.

Air-conditioning system recommended (room temperature inferior or equal to 25°C).

Material & Consumables necessary for the detection of genetic targets *via* real time PCR

ENRICHMENT :

Homogeneization

Stomacher

Stomacher bag with filter (precision : 400ml)

Incubation

Oven stabilised at 37°C or 30 °C

Transfer

Micropipettes : 25-200µl, 1 ml

Graded tube with cap (15ml)

Bacterial growth

Culture media (PCrone, EPT, Half Fraser ...)

Deionised water

Bacterial concentration (optional)

Graded tube with cap (1.5ml)

Centrifuge (rotor compatible for 1.5ml microtube)

DNA EXTRACTION

Pelleting reagents

Plate or strips centrifuge

Thermal and Chemical Lysis

Thermocycler or water bath or dry bath

Extraction buffer

Ethanol 96°

Transfer

Micropipette 25-200µl

PURIFICATION using COLUMNS

Pelleting reagents

Tubes, plates or strips centrifuge

Transfer

Micropipette 25-200µl and 300µl (1 or 8 channels)

Distripette with combitips

Filtration

Filtration Columns

4 Collecting tubes

Wash buffers 1 and 2

Elution

Elution buffer

Conic Tubes of 1.5ml volume

DNA Conservation

200µl Microtubes

PURIFICATION using MAGNETIC MICROBEADS

Pelleting reagents

Tubes, plates or strips centrifuge

Transfer

Micropipette 25-200µl and 300µl (1 or 8 channels)

Distripette with combitips

Reservoirs for buffers

Magnetic separation

Deep-well plate, X-pierced film, aluminium foil

Magnetic separating device 96

Wash buffers 1 and 2

Elution

Elution buffer

Conic Tubes of 1.5ml volume

DNA Conservation

AMPLIFICATION

Materials

Thermocycler (sybr channel)

Centrifuge with plate rotor

Plate holder

Microtube (200µl) plate holder

Reagents

Master mix PCR (1)

Master mix PCR Salmonella (AFNOR ADN 33/03 – 04/10) (1)

Master mix PCR O157 (AFNOR ADN 33/01 – 02/10) (1)

Master mix PCR H7 (AFNOR ADN 33/01 – 02/10) (1)

Master mix PCR Ti (inhibition control) (AFNOR ADN 33/01 – 02/10)

Standard DNA for each target (3µl)

(1) Strip of 8 colored wells filled beforehand, breakable, with individual optic caps.

2 Training

ADNucleis is registered for professional training :

SIREN number of the training organism: 497 802 165 000 17

TVA intraCEE : FR 3649780216500017

TRAINING activity registered under the
N° 82 69 09708 69 with Préfecture of Rhône-Alpes Region.

Fees and duration of training:

In accordance with the level of initial training of personnel and complexity of the demand.

Trainers schooling level : Master degree

Cost evaluation according to number of training days:

- For 1 fonctionnal laboratory : 1 work day
- For 1 laboratory in design: 5 work days

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